

**Isoflurane Consumption Chart**  
**Estimated CC's Per Hour**

Oxygen Flow (liters/min)	Vaporizer Setting					
	1.0	2.0	2.5	3.0	4.0	5.0
1.0	3.00	6.00	7.50	9.00	12.00	15.00
1.5	4.50	9.00	11.25	13.50	18.00	22.50
2.0	6.00	12.00	15.00	18.00	24.00	30.00
2.5	7.50	15.00	18.75	22.50	30.00	37.50
3.0	9.00	18.00	22.50	27.00	36.00	45.00

**Sevoflurane Consumption Chart**  
**Estimated CC's Per Hour**

Oxygen Flow (liters/min)	Vaporizer Setting					
	2.0	3.0	3.5	4.0	5.0	7.0
1.0	6.66	9.99	11.66	13.32	16.65	23.31
1.5	9.99	14.99	17.48	19.98	24.98	34.97
2.0	13.32	19.98	23.31	26.64	33.30	46.62
2.5	16.65	24.98	29.14	33.30	41.63	58.28
3.0	19.98	29.97	34.97	39.96	49.95	69.93